



# ElbaTech srl

Soluzioni elettroniche e Innovazione

## EQCM

Our top-quality QCMagic instrument can be now equipped with a precision Potentiostat directly integrated into the system.

Electrochemical experiments can thus be performed *concurrently* with QCM measurements.

This becomes possible because the Oscillator Unit exposes the quartz pad in contact with the measuring solution to the hard ground of the system.

For the potentiostat this pad thus behaves as the *Working Electrode*.

The user-friendly software offers both chronoamperometry measurements and cyclic voltammetry loops.

